


Double J Stent

U - 2026

- D. J. Stent used for temporary internal drainage from ureteropelvic junction to the bladder
- Made of superior grade poly-urethane material.
- Black line on the stent indicates the direction of coil after withdrawal of guide wire .
- Complete radio opaque stent for x-ray visualization. Implant tested with an externally low encrustation tendency.
- Ureteral Stent Consist : Stent, Pusher, with or without Thread.
- Ureteral Stent Kit Consist : Stent, Pusher, Guidewire, 2 Clamps with or without Thread.
Sterile / Disposable / Individually Packed.



Size in FG	4.0, 4.5, 4.8, 5.0, 5.5, 6.0
Length CM	16, 20, 26, 28
Type	Open end & Close end