

## Reinforced

### Cuffed & Plain

#### A - 3013

- Used to maintain a patent airway and mechanically ventilate patients.
  - Uniform spiral metal-reinforced tube wall to prevent the tube from kinking.
  - Reinforcement supports tube position and verification during X-ray.
  - Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line.
  - High volume / low pressure cuff ensures efficient low pressure cuff seal (For cuffed tubes).
  - Universal 15 mm connector at proximal end.
  - Available with DEHP free Material.
- Sterile / Disposable / Individually Packed.



	Size MM
CUFFED :	5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.00
PLAIN :	3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.00