



Bone Marrow Biopsy Needle

- Remover guide provided with indicators to facilitate the sample expulsion and to be able to check the length of the sample during the procedure.
- The "MM" and "MF" additional extraction cannulas allow to get the sample in a completely atraumatic and safe way and to hold it without any damage to its surface and inner structure until extraction.
- Available in Trocar, Aspiration or Fish Mouth tip.
- Luer Lock connection for syringe.

Sterile / Disposable / Individually Packed.

Color Code	Green	Pink	White	Light Blue
Size	8	9	11	13



