

Nebulizer Mask

A - 3015 - 3016

- Designed to be used with any standard Nebulizer kit.
 - Super smooth kink resistance tubing ensures uniform flowrate.
 - Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
 - Adjustable elastic strip and integrated nose clip for proper positioning of mask.
 - The Nebulizer Kit just plugs into the standard fitting located at the bottom of the Mask.
 - Masks are designed to deliver more medication directly to the patient due to the direct flow design.
 - Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
 - Star lumen main tube to avoid accidental blockage.
 - Tube Length : 210cm.
 - Size : Adult, Paediatric & Neonatal.
 - Available with DEHP free Material.
- Sterile / Disposable / Individually Packed.

