

T-PC. For Nebulizer

A - 3024



- Super smooth kink resistance tubing ensures uniform flowrate.
 - Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
 - Quick to assemble - Facilitates fast nebulisation
 - Can be used with Mask, iGel, LMA and ET tube
 - Rotatable sections - Can be used on patients in sitting or supine position.
 - Drug dosages and oxygen flow rates remain as per guidelines.
 - Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
 - Star lumen main tube to avoid accidental blockage.
 - Tube Length : 210cm.
 - Size : Adult, Paediatric & Neonatal.
 - Available with DEHP free Material.
- Sterile / Disposable / Individually Packed.